



Perennity DICOM edition
DICOM Conformance Statement

Version: 4.50

1. Table of contents

1.	Table of contents	2
2.	Preface	3
3.	Introduction	4
4.	Data workflow	5
5.	Perennity DICOM server (AE) specifications	6
5.1.	Transfer syntax	6
5.1.1.	Native	6
5.1.2.	Compression	6
5.1.3.	SR	6
5.2.	Supported SOP Classes	7
5.3.	Implementation Class and Version	11
6.	Verification	12
7.	C-STORE	13
8.	C-FIND	14
9.	C-MOVE	18
10.	C-ECHO	18

Revision history

Original document	Oct 2006	Laurent Janssens	
Revision 1	Oct 2009	Laurent Janssens	
Revision 2	May 2011	Laurent Janssens	
Rev.3 v4.10	May 2012	Laurent Janssens	
Rev. 4 v4.20	March 2013	Laurent Janssens	
Rev. 5 v4.21	Sep 2014	Laurent Janssens	
Rev. 5 v4.50	Aug 2015	Laurent Janssens	Review - Unchanged

2. Preface

Perennity DICOM edition offers a turnkey solution for the medical industry to output DICOM images from any modality or PACS system to CD, DVD and BluRay. While dramatically reducing production costs compared to traditional film, CD, DVD & BluRay can as well be used for archiving purposes as for the exchange of DICOM data (images, video, reports, ...) between hospitals, doctors and patients.

Perennity uses CD/DVD/BD autoloaders (robots) to fully automate the CD/DVD/BD recording and printing. Its operation can be 'attended', where a human action will trigger the writing processes, or 'unattended', in which case Perennity automatically generates the discs grouped by series, study or patient. IQ-Lite, Osirix, eFilm and Keosys DICOM viewers are available as an option and can be added to every CD, offering convenient browsing and image manipulation functionality (changing colors, windowing, measurements...). Those DICOM viewers, which autoruns when a CD is inserted, do not require any local installation on the PC and it can easily be translated in any language. Please note that Perennity also supports any other DICOM compatible viewers.

Perennity DICOM edition uses a central database (MySQL or MS SQL Server) to log all operations and keep track of which media contain which data. For archiving purposes you may optionally decide to keep available on the Perennity server all processed data - with the possibility to resubmit it to CD/DVD/BD later on.

Perennity DICOM is composed of two essential components; one being a DICOM Server module, the second one being the DICOM Burner module. Both are running as Windows services.

3. Introduction

This is a conformance statement for the "Perennity DICOM Server" and "Perennity DICOM Burner", which both support DICOM.


"Perennity DICOM Server" is a component of the Perennity DICOM edition solution. Its goal is to receive images sent by modalities or PACS system, and store them as file. It also extracts DICOM metadata in the Perennity database (MySQL or Microsoft SQL Server).

"Perennity DICOM Burner" is a second component of the Perennity DICOM solution. Its goal is to submit CD/DVD/BD Production jobs for images retrieved either from the local Perennity Archive or from remote PACS systems (via Query/Retrieve or Store from the PACS).

This document contains: SOP classes supported, "Perennity DICOM Server" roles and a view of interaction components in Perennity DICOM solution.

Supported roles: Storage SCP, Storage SCU, Query-Retrieve SCP

4. Data workflow



Different modalities and/or the PACS send images to the Perennity DICOM Server over the local area network using the DICOM Store protocol over TCP/IP, whereby the Perennity DICOM Server acts as a genuine DICOM server.

You may decide to keep all images in a high capacity storage (RAID, optical libraries, ...). If not, the images are temporarily stored in a buffer and will be deleted once processed. Images are automatically sorted by study and series. Perennity extracts their metadata (study ID, patient name or identifier, type of modality, etc.) and stores them in its database (MySQL or SQL Server).

Using the Perennity DICOM Burner application, the radiologist can query and select the series of images to be output to CD/DVD/BD. He can specify the number of copies and optionally add some comments to be printed on the disk and stored in the database.

The DICOM Burner can also be used to Query a third party PACS and Retrieve images from it.

Multiple studies belonging to different patients can also be combined on a single disk.

The CDs/DVDs/BDs are automatically recorded and printed using a CD/DVD/BD autoloader. This robot is hooked to a workstation or has a Windows based embedded PC. The robot is managed by the Perennity DICOM Burner.

The disc output process can also be fully automated, in which case every series will automatically be recorded to a CD/DVD/BD after a predefined latency time (e.g. when no images have been received during the last 5 minutes). A DICOM Worklist can also be used to select the studies/series to be recorded.

Series can be output to different discs or grouped by study or patient.

5. Perennity DICOM server (AE) specifications

5.1. Transfer syntax

5.1.1. Native

Name	UID
Explicit VR Little Endian	1.2.840.10008.1.2.1
Explicit VR Big Endian	1.2.840.10008.1.2.2
Implicit VR Little Endian	1.2.840.10008.1.2

5.1.2. Compression

Name	UID
JPEGBaseline	1.2.840.10008.1.2.4.50
JPEGExtended	1.2.840.10008.1.2.4.51
JPEGLSLossyNearLossless	1.2.840.10008.1.2.4.81
JPEG2000Lossy	1.2.840.10008.1.2.4.91
JPEGLossless14	1.2.840.10008.1.2.4.57
JPEGLossless	1.2.840.10008.1.2.4.70
JPEGLSLossless	1.2.840.10008.1.2.4.80
JPEG2000Lossless	1.2.840.10008.1.2.4.90
MPEG2	1.2.840.10008.1.2.4.100
RLELossless	1.2.840.10008.1.2.5

5.1.3. SR

Name	UID
Deflated Explicit VR Little Endian	1.2.840.10008.1.2.1.99

5.2. Supported SOP Classes

The Perennity DICOM server SCP entity provides standard conformance to the following SOP classes:

SOP Class Name	SOP Class UID	Transfer Syntax	Role
Basic Study Content Notification SOP Class (Retired)	1.2.840.10008.1.9	All	SCP
Stored Print Storage SOP Class (Retired)	1.2.840.10008.5.1.1.27	All	SCP
Hardcopy Grayscale Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.29	All	SCP
Hardcopy Color Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.30	All	SCP
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	All	SCP
Digital X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	All	SCP
Digital X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	All	SCP
Digital Mammography X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	All	SCP
Digital Mammography X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	All	SCP
Digital Intra-oral X-Ray Image Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.3	All	SCP
Digital Intra-oral X-Ray Image Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.3.1	All	SCP
Standalone Modality LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.10	All	SCP
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	All	SCP
Standalone VOI LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.11	All	SCP
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	All	SCP
Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.2	All	SCP
Pseudo-Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.3	All	SCP

SOP Class Name	SOP Class UID	Transfer Syntax	Role
Blending Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.4	All	SCP
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	All	SCP
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1	All	SCP
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	All	SCP
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1	All	SCP
X-Ray Angiographic Bi-Plane Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.12.3	All	SCP
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	All	SCP
Standalone PET Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.129	All	SCP
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	All	SCP
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1	All	SCP
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	All	SCP
Ultrasound Multi-frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3	All	SCP
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	All	SCP
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	All	SCP
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1	All	SCP
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2	All	SCP
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	All	SCP
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	All	SCP
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	All	SCP
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4	All	SCP
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	All	SCP
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6	All	SCP
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7	All	SCP

SOP Class Name	SOP Class UID	Transfer Syntax	Role
Nuclear Medicine Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5	All	SCP
Ultrasound Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6	All	SCP
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	All	SCP
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66	All	SCP
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1	All	SCP
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2	All	SCP
Real World Value Mapping Storage	1.2.840.10008.5.1.4.1.1.67	All	SCP
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	All	SCP
Multi-frame Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1	All	SCP
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	All	SCP
Multi-frame Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3	All	SCP
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	All	SCP
VL Image Storage - Trial (Retired)	1.2.840.10008.5.1.4.1.1.77.1	All	SCP
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	All	SCP
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1	All	SCP
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	All	SCP
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1	All	SCP
VL Slide-Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	All	SCP
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	All	SCP
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1	All	SCP

SOP Class Name	SOP Class UID	Transfer Syntax	Role
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	All	SCP
Ophthalmic Photography 16 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2	All	SCP
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3	All	SCP
VL Multi-frame Image Storage - Trial (Retired)	1.2.840.10008.5.1.4.1.1.77.2	All	SCP
Standalone Overlay Storage (Retired)	1.2.840.10008.5.1.4.1.1.8	All	SCP
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	All	SCP
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	All	SCP
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	All	SCP
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40	All	SCP
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	All	SCP
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	All	SCP
Chest CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.65	All	SCP
Standalone Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.9	All	SCP
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	All	SCP
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	All	SCP
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	All	SCP
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	All	SCP
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	All	SCP
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	All	SCP
Hanging Protocol Storage	1.2.840.10008.5.1.4.38.1	All	SCP
Siemens CSA Non-Image Storage	1.3.12.2.1107.5.9.1	All	SCP

5.3. Implementation Class and Version

Implementation Class UID	1.2.40.0.13.1.1
Implementation Version Name	Dcm4che-2.0

6. Verification

SOP Class Name	SOP Class UID	Role
Verification	1.2.840.10008.1.1	SCP

7. C-STORE

SOP Class Name	SOP Class UID	Role
Storage	1.2.840.10008.4.2	SCP

“Perennity DICOM Server” implements C-STORE service like defined in standard DICOM (ASSOCIATION-RQ, A-ASSOCIATE-AC, etc). The SOP Class, which can be accepted for the transfer, are included in the “Supported SOP classes” table of this document.

Service status	Return code
Success	0000
Error	0001

8. C-FIND

SOP Class Name	SOP Class UID	Role
Patient Root Query/Retrieve information Model – FIND	1.2.840.10008.5.1.4.1.2.1.1	SCP
Study Root Query/Retrieve information Model – FIND	1.2.840.10008.5.1.4.1.2.2.1	SCP
Patient/Study Only Query/Retrieve information Model – FIND	1.2.840.10008.5.1.4.1.2.3.1	SCP

Function C-FIND allows querying the “Perennity DICOM Server” to obtain information related to the images already stored. Four retrieve levels are available: PATIENT, STUDY, SERIES or IMAGES.

Information returned for each item:

PATIENT

Attribute Name	Tag	VR	Request	Return
<i>SOP Common</i>				
Specific Character Set	0008,0005	CS	x	x
<i>Patient Level</i>				
Patient's Name	0010,0010	PN	x	x
Patient ID	0010,0020	LO	x	x
Patient's Birth Date	0010,0030	DA	x	x
Patient's Sex	0010,0040	CS	x	x
Number Of Patient Related Studies	0020,1200	IS		x
Number Of Patient Related Series	0020,1202	IS		x
Number Of Patient Related Instances	0020,1204	IS		x

STUDY

Attribute Name	Tag	VR	Request	Return
<i>SOP Common</i>				
Specific Character Set	0008,0005	CS	x	x
<i>Patient Level</i>				
Patient's Name	0010,0010	PN	x	x
Patient ID	0010,0020	LO	x	x
Patient's Birth Date	0010,0030	DA	x	x
Patient's Sex	0010,0040	CS	x	x
<i>Study Level</i>				
Study Instance UID	0020,000D	UI	x	x
Study ID	0020,0010	SH	x	x
Study Date	0008,0020	DA	x	x
Study Description	0008,1030	LO	x	x
Accession Number	0008,0050	SH	x	x
Referring Physician's Name	0008,0090	PN	x	x
Modalities In Study	0008,0061	CS	x	x
Institution Name	0008,0080	LO		x
Number Of Study Related Series	0020,1206	IS		x
Number Of Study Related Instances	0020,1208	IS		x

SERIES

Attribute Name	Tag	VR	Request	Return
<i>SOP Common</i>				
Specific Character Set	0008,0005	CS	x	x
<i>Patient Level</i>				
Patient's Name	0010,0010	PN	x	
Patient ID	0010,0020	LO	x	
Patient's Birth Date	0010,0030	DA	x	
Patient's Sex	0010,0040	CS	x	
<i>Study Level</i>				
Study Instance UID	0020,000D	UI	x	
Study ID	0020,0010	SH	x	
Study Date	0008,0020	DA	x	
Study Description	0008,1030	LO	x	
Accession Number	0008,0050	SH	x	
Referring Physician's Name	0008,0090	PN	x	
Modalities In Study	0008,0061	CS	x	
Modality	0008,0060	CS	x	x
Series Number	0020,0011	IS	x	x
Series Instance UID	0020,000E	UI	x	x
Series Description	0008,103E	LO		x
Series Date	0008,0021	DA		x
Number Of Series Related Instances	0020,1209	IS		x

SERIES

Attribute Name	Tag	VR	Request	Return
<i>SOP Common</i>				
Specific Character Set	0008,0005	CS	x	x
<i>Patient Level</i>				
Patient's Name	0010,0010	PN	x	
Patient ID	0010,0020	LO	x	
Patient's Birth Date	0010,0030	DA	x	
Patient's Sex	0010,0040	CS	x	
<i>Study Level</i>				
Study Instance UID	0020,000D	UI	x	
Study ID	0020,0010	SH	x	
Study Date	0008,0020	DA	x	
Study Description	0008,1030	LO	x	
Accession Number	0008,0050	SH	x	
Referring Physician's Name	0008,0090	PN	x	
Modalities In Study	0008,0061	CS	x	
Modality	0008,0060	CS	x	
Series Number	0020,0011	IS	x	
Series Instance UID	0020,000E	UI	x	
Instance Number	0020,0013	IS	x	x
SOP Class UID	0008,0016	UI	x	x
SOP Instance UID	0008,0018	UI	x	x

The complete set of data is returned during the C-FIND query. If the SCU specifies a date (or a range of two dates – from ... to ...) in the "StudyDate" tag, the result set will be refined and returns only the data within this date range.

9. C-MOVE

SOP Class Name	SOP Class UID	Role
Patient Root Query/Retrieve information Model – MOVE	1.2.840.10008.5.1.4.1.2.1.2	SCP
Study Root Query/Retrieve information Model – MOVE	1.2.840.10008.5.1.4.1.2.2.2	SCP
Patient/Study Only Query/Retrieve information Model – MOVE	1.2.840.10008.5.1.4.1.2.3.2	SCP

10. C-ECHO

Service Status	Return Code	Behavior
Success	0000	The C-ECHO request is accepted